

Anti-Icam-1/Cd54 Mouse Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P13597

Item Number: CM0008850

Introduction

High-specificity anti-mouse ICAM-1/CD54 mouse monoclonal antibody validated for WB, IF-P, IHC-P, and ELISA. Recognizes mouse and rat ICAM-1 for cell adhesion and leukocyte migration studies.

[Learn More](#)

Parameter	Value
Product Name	ICAM-1/CD54 Mouse mAb
Remarks/Alias	CD54; Ly-47; Icam-1; MALA-2
Species	Mouse
GeneID (Human)	3383
GeneID	15894
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 28-485 of mouse ICAM-1/CD54 (NP_034623.1)
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P13597
Protein Weight	59 kDa
Shipping	Ice bag