


Anti-Icam-1/Cd54 Polyclonal Antibody For Wb And Elisa - P05362

Item Number: CM0017515



Introduction

High-quality rabbit polyclonal antibody against human ICAM-1/CD54 (BB2, CD54), validated for WB and ELISA. Cross-reacts with human and mouse. Suitable for adhesion, inflammation, and viral entry research. UniProt: P05362.

[Learn More](#)

Parameter	Value
Product Name	ICAM-1/CD54 Rabbit pAb
Remarks/Alias	BB2; CD54; P3.58; ICAM-1/CD54
Species	Human
GeneID (Human)	3383
GeneID	3383
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 250-480 of human ICAM-1/CD54 (NP_000192.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P05362
Protein Weight	58kDa
Shipping	Ice bag