

Anti-Icam-1/Cd54 Rabbit Monoclonal Antibody For Wb, Ihc-P, Fc, Elisa - Q00238

Item Number: CM0011796

Introduction

Rabbit monoclonal antibody against rat ICAM-1/CD54, validated for WB, IHC-P, FC, ELISA. Reacts with rat. Ideal for studying leukocyte adhesion, trans-endothelial migration, and inflammation.

[Learn More](#)

Parameter	Value
Product Name	ICAM-1/CD54 Rabbit mAb
Remarks/Alias	CD54; ICAM
Species	Rat
GenelD (human)	3383
GenelD	25464
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-493 of rat ICAM-1/CD54 (NP_037099.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, FC, ELISA
Cross Reactivity	Rat
SWISS	Q00238
Protein Weight	60kDa
Shipping	Ice bag