


Anti-Icam-1/Cd54 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P05362

Item Number: CM0027566



Introduction

Rabbit polyclonal antibody targeting human ICAM-1/CD54 (SWISS P05362). Validated for WB, IHC-P, IF/ICC, ELISA with cross-reactivity to mouse and rat. Ideal for studying cell adhesion, leukocyte migration, and viral infection mechanisms.

[Learn More](#)

Parameter	Value
Product Name	ICAM-1/CD54 Rabbit pAb
Remarks/Alias	BB2; CD54; P3.58; ICAM-1/CD54
Species	Human
GeneID (Human)	3383
GeneID	3383
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-480 of human ICAM-1/CD54 (NP_000192.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05362
Protein Weight	58kDa
Shipping	Ice bag