

Anti-Icam1 (Cd54) Polyclonal Antibody For Wb, If-P, Elisa - Q00238

Item Number: CM0023438

Introduction

Icam1 Rabbit pAb (polyclonal antibody) validated for WB, IF-P, ELISA with cross-reactivity to mouse and rat. Targets intercellular adhesion molecule 1 (CD54), 60kDa, UniProt Q00238. Ideal for cell adhesion and leukocyte migration studies.

[Learn More](#)

Parameter	Value
Product Name	Icam1 Rabbit pAb
Remarks/Alias	CD54; ICAM; Icam1
Species	Rat
GeneID (Human)	3383
GeneID	25464
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-493 of rat Icam1 (NP_037099.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q00238
Protein Weight	60kDa
Shipping	Ice bag