

Anti-Icam1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P13597

Item Number: CM0021886

Introduction

High-quality rabbit polyclonal antibody against mouse ICAM1 (CD54), validated for WB, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for leukocyte adhesion and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	ICAM1 Rabbit pAb
Remarks/Alias	CD54; Ly-47; Icam-1; MALA-2; ICAM1
Species	Mouse
GeneID (human)	3383
GeneID	15894
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 120-320 of mouse ICAM1 (NP_034623.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13597
Protein Weight	59kDa
Shipping	Ice bag