

Anti-Icam2 / Cd102 Rabbit Monoclonal Antibody For Wb, Elisa - P35330

Item Number: CM0003157



Introduction

Recombinant anti-mouse ICAM2/CD102 rabbit monoclonal antibody for Western blot and ELISA. Cross-reacts with mouse and rat. Ideal for studying leukocyte adhesion and immune cell interactions.

[Learn More](#)

Parameter	Value
Product Name	ICAM2 / CD102 Rabbit mAb
Remarks/Alias	CD102; Ly-60; Icam-2; ICAM2 / CD102
Species	Mouse
GeneID	15896
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of mouse ICAM2 / CD102 (P35330)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P35330
Protein Weight	31kDa
Shipping	Ice bag