

Anti-Icam2/Cd102 Rabbit Polymab® For Wb, Ihc-P, Elisa - P35330

Item Number: CM0006114



Introduction

Anti-ICAM2/CD102 Rabbit PolymAb® (recombinant monoclonal) validated for WB, IHC-P, ELISA. Recognizes mouse ICAM2, a cell adhesion molecule involved in leukocyte migration. Recombinant immunogen, 31kDa, shipping on ice.

[Learn More](#)

Parameter	Value
Product Name	ICAM2/CD102 Rabbit PolymAb®
Remarks/Alias	CD102; Ly-60; Icam-2
Species	Mouse
GeneID (Human)	3384
GeneID	15896
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-222 of mouse ICAM2/CD102 (NP_034624.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse
SWISS	P35330
Protein Weight	31kDa
Shipping	Ice bag