

Anti-Icam3 Rabbit Polyclonal Antibody For Wb, Elisa - P32942

Item Number: CM0024941



Introduction

Rabbit polyclonal antibody against human ICAM3 (CD50), validated for WB and ELISA. Ideal for studying leukocyte adhesion, phagocytosis, and dendritic cell interactions. Cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	ICAM3 Rabbit pAb
Remarks/Alias	CD50; CDW50; ICAM-R; ICAM3
Species	Human
Gene ID (Human)	3385
Gene ID	3385
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 270-490 of human ICAM3 (NP_002153.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P32942
Protein Weight	60kDa
Shipping	Ice bag