

Anti-Icam3/Cd50 Rabbit Monoclonal Antibody For Wb, Ihc-P, Fc - P32942

Item Number: CM0009460



Introduction

Rabbit monoclonal antibody against human ICAM3/CD50. Validated for WB, IF/ICC, IHC-P, FC, ELISA. Cross-reacts with human. Ideal for studying leukocyte adhesion and immune cell interactions.

[Learn More](#)

Parameter	Value
Product Name	ICAM3/CD50 Rabbit mAb
Remarks/Alias	CD50; CDW50; ICAM-R
Species	Human
GeneID (Human)	3385
GeneID	3385
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-114 of human ICAM3 (NP_002153.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, FC, ELISA
Cross Reactivity	Human
SWISS	P32942
Protein Weight	60kDa
Shipping	Ice bag