

Anti-Icam4 Polyclonal Antibody For Wb, Elisa - Q14773

Item Number: CM0029058

Introduction

High-quality rabbit polyclonal antibody against ICAM4 (LW blood group protein). Validated for Western blot and ELISA, cross-reacts with human, mouse, and rat. Ideal for erythrophagocytosis and platelet adhesion studies.

[Learn More](#)

Parameter	Value
Product Name	ICAM4 Rabbit pAb
Remarks/Alias	LW; CD242; ICAM4
Species	Human
GeneID (Human)	3386
GeneID	3386
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-140 of human ICAM4 (NP_001034221.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14773
Protein Weight	29kDa
Shipping	Ice bag