

## Anti-Ick Polyclonal Antibody For Wb, Elisa - Q9Upz9

Item Number: CM0026516

### Introduction

IVD-ready anti-ICK polyclonal antibody validated for WB and ELISA. Cross-reacts with human, mouse, and rat. Targets serine/threonine-protein kinase ICK (Q9UPZ9) for ciliogenesis, developmental, and cancer research.

[Learn More](#)

| Parameter        | Value                                                                                                                               |
|------------------|-------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | ICK Rabbit pAb                                                                                                                      |
| Remarks/Alias    | ECO; ICK; MRK; LCK2; hICK; EJM10                                                                                                    |
| Species          | Human                                                                                                                               |
| Gene ID (Human)  | 22858                                                                                                                               |
| Gene ID          | 22858                                                                                                                               |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human ICK (NP_057597.2). |
| Source           | Rabbit                                                                                                                              |
| Category         | Polyclonal Antibodies                                                                                                               |
| Application      | WB, ELISA                                                                                                                           |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                   |
| SWISS            | Q9UPZ9                                                                                                                              |
| Protein Weight   | 71kDa                                                                                                                               |
| Shipping         | Ice bag                                                                                                                             |