

Anti-Icmt Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O60725

Item Number: CM0015769

Introduction

High-quality ICMT rabbit polyclonal antibody validated for WB, IHC-P, and ELISA in human, mouse, and rat. Ideal for studying protein prenylation and methylation pathways.

[Learn More](#)

Parameter	Value
Product Name	ICMT Rabbit pAb
Remarks/Alias	PCMT; PPMT; PCCMT; HSTE14; MST098; MSTP098; ICMT
Species	Human
GeneID (Human)	23463
GeneID	23463
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 175-284 of human ICMT (NP_036537.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60725
Protein Weight	32kDa
Shipping	Ice bag