


Anti-Icos/Cd278 Monoclonal Antibody For Wb And Elisa - Q9Y6W8

Item Number: CM0012342



Introduction

Rabbit monoclonal antibody against human ICOS/CD278, validated for WB and ELISA. Cross-reacts with mouse and rat. Targets a key T-cell costimulatory receptor critical for immune responses and antibody production.

[Learn More](#)

Parameter	Value
Product Name	ICOS/CD278 Rabbit mAb
Remarks/Alias	AILIM; CD278; CVID1; ICOS/CD278
Species	Human
GeneID	29851
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ICOS/CD278 (Q9Y6W8)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y6W8
Protein Weight	23kDa
Shipping	Ice bag