

Anti-Icos/Cd278 Rabbit Monoclonal Antibody For Wb, If/Icc, Fc, Elisa - Q9Wvs0

Item Number: CM0006708

Introduction

Rabbit monoclonal antibody against mouse ICOS (CD278), validated for WB, IF/ICC, FC, and ELISA. Suitable for T-cell costimulation and immune response research. Reacts with mouse. SWISS: Q9WVS0.

[Learn More](#)

Parameter	Value
Product Name	ICOS/CD278 Rabbit mAb
Remarks/Alias	H4; CCLP; AILIM; CRP-1; Ly115; ICOS/CD278
Species	Mouse
GeneID (Human)	29851
GeneID	54167
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-142 of mouse ICOS/CD278(NP_059508.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, FC, ELISA
Cross Reactivity	Mouse
SWISS	Q9WVS0
Protein Weight	22kDa
Shipping	Ice bag