

Anti-Icosl Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O75144

Item Number: CM0026343



Introduction

High-quality rabbit polyclonal antibody against human ICOSL (CD275/B7-H2). Validated for Western blot, IHC-P, and ELISA. Cross-reacts with mouse and rat. Suitable for T-cell costimulation and immune checkpoint studies.

[Learn More](#)

Parameter	Value
Product Name	ICOSL Rabbit pAb
Remarks/Alias	B7h, B7H2, GL50, B7-H2, B7RP1, CD275, ICOSL, LICOS, B7RP-1, ICOS-L
Species	Human
GeneID (Human)	23308
GeneID	23308
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-256 of human ICOSL (NP_056074.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75144
Protein Weight	33kDa
Shipping	Ice bag