


Anti-Ict1/Mrp-L58 Rabbit Polyclonal Antibody For Wb, Elisa - Q14197

Item Number: CM0016969



Introduction

Anti-ICT1/MRP-L58 rabbit polyclonal antibody for WB and ELISA detection of human and mouse ICT1, a mitochondrial peptidyl-tRNA hydrolase. Ideal for mitochondrial translation termination and ribosome recycling research. – Q14197

[Learn More](#)

Parameter	Value
Product Name	ICT1 Rabbit pAb
Remarks/Alias	DS1; DS-1; ICT1; MRP-L58
Species	Human
Gene ID (Human)	3396
Gene ID	3396
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 30-206 of human ICT1 (NP_001536.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q14197
Protein Weight	24kDa
Shipping	Ice bag