

Anti-Id1 Rabbit Polyclonal Antibody For Wb & Elisa - P41134

Item Number: CM0026351

Introduction

High-quality anti-ID1 rabbit polyclonal antibody for WB and ELISA applications. Validated in human, mouse, and rat samples. Targets DNA-binding protein inhibitor ID-1 with a molecular weight of 16 kDa. Ideal for studying transcriptional regulation, cell differentiation, and circadian rhythms.

[Learn More](#)

Parameter	Value
Product Name	ID1 Rabbit pAb
Remarks/Alias	ID; bHLHb24; ID1
Species	Human
GeneID (human)	3397
GeneID	3397
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-155 of human ID1 (NP_002156.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P41134
Protein Weight	16kDa
Shipping	Ice bag