

Anti-Id2 Rabbit Polyclonal Antibody For Wb, Elisa - Q02363

Item Number: CM0026716



Introduction

Rabbit polyclonal antibody against human ID2, validated for WB and ELISA. Targets the transcriptional regulator ID2, a bHLH inhibitor involved in differentiation, circadian rhythm, and cancer. UniProt: Q02363.

[Learn More](#)

Parameter	Value
Product Name	ID2 Rabbit pAb
Remarks/Alias	GIG8; ID2A; ID2H; bHLHb26; ID2
Species	Human
GeneID (Human)	3398
GeneID	3398
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ID2 (NP_002157.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q02363
Protein Weight	15kDa
Shipping	Ice bag