

Anti-Id3 Rabbit Polyclonal Antibody For Ihc-P, If/Icc, Elisa - Q02535

Item Number: CM0027711



Introduction

Rabbit polyclonal antibody targeting human ID3, suitable for IHC-P, IF/ICC, ELISA. Recognizes human, mouse, rat ID3. Investigate transcriptional regulation, myogenesis, and circadian clock pathways.

[Learn More](#)

Parameter	Value
Product Name	ID3 Rabbit pAb
Remarks/Alias	HEIR-1; bHLHb25; ID3
Species	Human
Gene ID (Human)	3399
Gene ID	3399
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-119 of human ID3 (NP_002158.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02535
Protein Weight	13kDa
Shipping	Ice bag