

Anti-Ide Rabbit Monoclonal Antibody [Kd Validated] For Wb, Ihc-P, Elisa - P14735

Item Number: CM0004578



Introduction

High-quality rabbit monoclonal antibody against human IDE, validated for WB, IHC-P, and ELISA. Recognizes human, mouse, and rat IDE. Ideal for studying insulin degradation and amyloid-beta clearance. SWISS: P14735.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] Insulin-degrading enzyme/IDE Rabbit mAb
Remarks/Alias	INSULYSIN
Species	Human
Gene ID (Human)	3416
Gene ID	3416
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human Insulin-degrading enzyme (IDE) (NP_004960.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14735
Protein Weight	118kDa
Shipping	Ice bag