

Anti-Idh1 Polyclonal Antibody For Wb, If/Icc, Elisa - 075874

Item Number: CM0024338



Introduction

CamelBio's IDH1 Rabbit pAb targets cytoplasmic isocitrate dehydrogenase, validated for WB, IF/ICC, ELISA in human, mouse, rat, monkey. Essential for studying cancer metabolism and TCA cycle.

[Learn More](#)

Parameter	Value
Product Name	IDH1 Rabbit pAb
Remarks/Alias	IDH; IDP; IDCD; IDPC; PICD; HEL-216; HEL-S-26; IDH1
Species	Human
GeneID (Human)	3417
GeneID	3417
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-200 of human IDH1 (NP_005887.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Monkey
SWISS	075874
Protein Weight	47kDa
Shipping	Ice bag