

Anti-Idh1 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - 075874

Item Number: CM0007639



Introduction

Rabbit monoclonal antibody targeting human IDH1, validated for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Detects cytoplasmic NADP-isocitrate dehydrogenase critical for NADPH generation and metabolic pathways.

[Learn More](#)

Parameter	Value
Product Name	IDH1 Rabbit mAb
Remarks/Alias	IDH; IDP; IDCD; IDPC; PICD; HEL-216; HEL-S-26; IDH1
Species	Human
Gene ID (Human)	3417
Gene ID	3417
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075874
Protein Weight	47kDa
Shipping	Ice bag