

Anti-Idh2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa, Chip - P48735

Item Number: CM0024866



Introduction

Rabbit polyclonal antibody targeting human IDH2, a mitochondrial TCA cycle enzyme. Validated in WB, IHC-P, IF/ICC, IP, ELISA, ChIP; cross-reacts with mouse and rat. Ideal for metabolism and cancer studies.

[Learn More](#)

Parameter	Value
Product Name	IDH2 Rabbit pAb
Remarks/Alias	IDH; IDP; IDHM; IDPM; ICD-M; IDH-2; D2HGA2; mNADP-IDH; IDH2
Species	Human
GeneID (human)	3418
GeneID	3418
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 193-452 of human IDH2 (NP_002159.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat
SWISS	P48735
Protein Weight	51kDa
Shipping	Ice bag