

Anti-Idh3B Rabbit Polyclonal Antibody For Wb, If, Elisa - O43837

Item Number: CM0017369



Introduction

Rabbit polyclonal antibody against IDH3B, the mitochondrial NAD-dependent isocitrate dehydrogenase beta subunit. Validated for WB, IF/ICC, and ELISA with human, mouse, and rat reactivity. Ideal for TCA cycle and mitochondrial function research.

[Learn More](#)

Parameter	Description
Product Name	IDH3B Rabbit pAb
Remarks/Alias	RP46; IDH3B
Species	Human
Gene ID (Human)	3420
Gene ID	3420
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-170 of human IDH3B (NP_008830.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43837
Protein Weight	42kDa
Shipping	Ice bag