


Anti-Ido1 Polyclonal Antibody For Wb, If/Icc, Elisa - P14902

Item Number: CM0019037



Introduction

Rabbit polyclonal antibody targeting human IDO1, validated for Western blot, immunofluorescence, and ELISA. Detects human and rat IDO1. Essential tool for tryptophan metabolism and immune regulation research.

[Learn More](#)

Parameter	Value
Product Name	IDO1 Rabbit pAb
Remarks/Alias	IDO; INDO; IDO-1; IDO1
Species	Human
Gene ID	3620 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-403 of human IDO1 (NP_002155.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	P14902
Protein Weight	45kDa
Shipping	Ice bag