


Anti-Ido1 Rabbit Polymab® For Wb, If/Icc, Elisa - P14902

Item Number: CM0003468



Introduction

High-quality Anti-IDO1 Rabbit PolymAb® for WB, IF/ICC, ELISA applications. Detects human IDO1 (P14902), a key enzyme in tryptophan catabolism and immune tolerance, suitable for cancer and immunology research.

[Learn More](#)

Parameter	Value
Product Name	IDO1 Rabbit PolymAb®
Remarks/Alias	IDO; INDO; IDO-1
Species	Human
GeneID (Human)	3620
GeneID	3620
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-403 of human IDO1 (NP_002155.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P14902
Protein Weight	45kDa
Shipping	Ice bag