


Anti-Ids Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P22304

Item Number: CM0023420



Introduction

Rabbit polyclonal antibody targeting human Iduronate 2-sulfatase (IDS). Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for lysosomal storage disease research and biomarker studies.

[Learn More](#)

Specification	Details
Product Name	IDS Rabbit pAb
Remarks/Alias	ID2S; MPS2; SIDS; IDS
Species	Human
GeneID (Human)	3423
GeneID	3423
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 34-290 of human IDS (NP_000193.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22304
Protein Weight	62kDa
Shipping	Ice bag