

Anti-Ier5 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q5Vy09

Item Number: CM0020744



Introduction

Rabbit polyclonal antibody against immediate early response gene 5 protein (IER5), validated for WB, IHC-P, IF/ICC, ELISA. Suitable for human, mouse, rat samples. Involved in transcription regulation, heat shock stress response, and cell cycle control.

[Learn More](#)

Parameter	Value
Product Name	IER5 Rabbit pAb
Remarks/Alias	SBBI48; IER5
Species	Human
GeneID (Human)	51278
GeneID	51278
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 200-300 of human IER5 (NP_057629.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5VY09
Protein Weight	34kDa
Shipping	Ice bag