

Anti-IFI16 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q16666

Item Number: CM0017380



Introduction

Rabbit polyclonal antibody against human IFI16 (alias PYHIN2, IFNGIP1) validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Useful for innate immunity, inflammation, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	IFI16 Rabbit pAb
Remarks/Alias	PYHIN2; IFNGIP1; IFI16
Species	Human
Gene ID (Human)	3428
Gene ID	3428
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human IFI16 (NP_001351796.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16666
Protein Weight	88kDa
Shipping	Ice bag