

Anti-Ifi16 Rabbit Polyclonal Antibody - Q16666

Item Number: CM0022818



Introduction

IFI16 rabbit polyclonal antibody for WB, FC, ELISA. Detects human IFI16, an innate immunity DNA sensor and transcriptional regulator involved in viral defense and cancer. Ideal for immunodetection studies.

[Learn More](#)

Parameter	Value
Product Name	IFI16 Rabbit pAb
Remarks/Alias	IFI16; IFNGIP1; PYHIN2; gamma-interferon-inducible protein 16
Species	Human
Gene ID (Human)	3428
Gene ID	3428
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 500-785 of human IFI16 (NP_001351796.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Human
SWISS	Q16666
Protein Weight	75kDa/82kDa/88kDa
Shipping	Ice bag