

Anti-Ifi27 Rabbit Polyclonal Antibody For Wb And Elisa - P40305

Item Number: CM0015947

Introduction

The IFI27 rabbit polyclonal antibody (alias ISG12A) specifically recognizes interferon alpha-inducible protein 27, mitochondrial. Validated for Western blot and ELISA applications, it cross-reacts with rat. Ideal for investigating type-I interferon-induced apoptosis, viral defense mechanisms, and mitochondrial signaling. Generated against a synthetic peptide within amino acids 1-100 of human IFI27.

[Learn More](#)

Parameter	Value
Product Name	IFI27 Rabbit pAb
Remarks/Alias	P27; ISG12; FAM14D; ISG12A; IFI27
Species	Human
GeneID	3429 (Human)
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human IFI27 (NP_005523.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P40305
Protein Weight	12kDa
Shipping	Ice bag