

Anti-Ifi30 Rabbit Polyclonal Antibody For Wb And Elisa - P13284

Item Number: CM0015470

Introduction

Rabbit polyclonal antibody against human IFI30 (GILT), suitable for Western blot and ELISA. Useful for studying MHC class I/II antigen processing and lysosomal thiol reductase activity in antigen-presenting cells.

[Learn More](#)

Parameter	Value
Product Name	IFI30 Rabbit pAb
Remarks/Alias	GILT; IP30; IP-30; IFI-30; IFI30
Species	Human
GeneID (Human)	10437
GeneID	10437
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human IFI30 (NP_006323.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P13284
Protein Weight	28kDa
Shipping	Ice bag