

Anti-Ifi35 Rabbit Polyclonal Antibody For Wb And Elisa - P80217

Item Number: CM0018777

Introduction

Rabbit polyclonal antibody against human IFI35 (IFP35), validated for Western blot and ELISA. Reacts with human and mouse. IFI35 is an interferon-induced protein regulating innate immune responses and NF-kappa-B signaling.

[Learn More](#)

Parameter	Value
Product Name	IFI35 Rabbit pAb
Remarks/Alias	IFP35; IFI35
Species	Human
GeneID (Human)	3430
GeneID	3430
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-288 of human IFI35 (NP_005524.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P80217
Protein Weight	32kDa
Shipping	Ice bag