

Anti-Ifi44 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q8Tcb0

Item Number: CM0027145

Introduction

Rabbit polyclonal antibody targeting IFI44, validated for WB, IHC-P, IF/ICC, and ELISA. Recognizes human, mouse, and rat IFI44. Suitable for interferon-induced microtubular structure studies and proteomics research.

[Learn More](#)

Parameter	Value
Product Name	IFI44 Rabbit pAb
Remarks/Alias	p44; MTAP44; A430056A10Rik; IFI44
Species	Mouse
GeneID (Human)	10561
GeneID	99899
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of mouse IFI44 (NP_598632.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TCB0
Protein Weight	50kDa
Shipping	Ice bag