


Anti-Ifit2 Polyclonal Antibody For Wb, Elisa - P09913

Item Number: CM0012989



Introduction

Rabbit polyclonal antibody targeting human IFIT2 (ISG54), validated for WB and ELISA. Suitable for antiviral immunity research, interferon signaling, and apoptosis studies.

[Learn More](#)

Parameter	Details
Product Name	IFIT2 Rabbit pAb
Remarks/Alias	P54, G10P2, IFI54, ISG54, cig42, IFI-54, IFIT-2, ISG-54, GARG-39, IFI-54K, ISG-54K, ISG-54 K, IFIT2
Species	Human
Gene ID (Human)	3433
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 273-472 of human IFIT2 (NP_001538.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P09913
Protein Weight	55kDa
Shipping	Ice bag