

Anti-Ifit3 Polyclonal Antibody For Wb, Elisa - O14879

Item Number: CM0015713



Introduction

Anti-IFIT3 polyclonal antibody validated for WB and ELISA, suitable for human, mouse, and rat samples. Explore antiviral and innate immune research with this high-quality reagent from CamelBio.

[Learn More](#)

Parameter	Value
Product Name	IFIT3 Rabbit pAb
Remarks/Alias	P49; Ifi49; IFIT3
Species	Mouse
GeneID (Human)	3437
GeneID	15959
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of mouse IFIT3 (NP_034631.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14879
Protein Weight	56kDa
Shipping	Ice bag