

Anti-Ifit3 Polyclonal Antibody For Wb, Ihc-P, Elisa - O14879

Item Number: CM0030126



Introduction

IFIT3 rabbit polyclonal antibody validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat IFIT3. Suitable for studying innate antiviral immunity, cell cycle regulation, and interferon signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	IFIT3 Rabbit pAb
Remarks/Alias	P49, Ifi49, IFIT3
Species	Mouse
GeneID (Human)	3437
GeneID	15959
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 128-331 of mouse IFIT3 (NP_034631.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14879
Protein Weight	56kDa
Shipping	Ice bag