

Anti-Ifitm3 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q01628

Item Number: CM0006598

Introduction

High-quality rabbit monoclonal antibody targeting human IFITM3 (1-8U, IP15, DSPA2b), validated for WB, IHC-P, and ELISA. Ideal for studying antiviral innate immunity and viral entry mechanisms.

[Learn More](#)

Parameter	Value
Product Name	IFITM3 Rabbit mAb
Remarks/Alias	1-8U; IP15; DSPA2b; IFITM3
Species	Human
GeneID (human)	10410
GeneID	10410
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-57 of human IFITM3(NP_066362.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q01628
Protein Weight	15kDa
Shipping	Ice bag