


Anti-Ifn-Beta Polyclonal Antibody For Wb, If/Icc, If-P, Elisa - P01575

Item Number: CM0020872



Introduction

Rabbit polyclonal antibody targeting mouse IFN-beta (UniProt P01575), validated for WB, IF/ICC, IF-P, and ELISA. Cross-reacts with human and mouse. Ideal for innate immunity and interferon signaling research.

[Learn More](#)

Parameter	Value
Product Name	IFN-beta Rabbit pAb
Remarks/Alias	Ifb; IFNB; If1da1; IFN-beta
Species	Mouse
GeneID (Human)	3456
GeneID	15977
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-182 of mouse IFN-beta (NP_034640.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P01575
Protein Weight	22kDa
Shipping	Ice bag