


Anti-Ifn-Beta Rabbit Monoclonal Antibody For Wb, Elisa - P01574

Item Number: CM0012056



Introduction

Rabbit monoclonal antibody against human IFN-beta for Western blot and ELISA applications. Suitable for innate immunity, antiviral response, and cancer research. Recombinant immunogen, detects 22 kDa protein.

[Learn More](#)

Parameter	Value
Product Name	IFN-beta Rabbit mAb
Remarks/Alias	IFB; IFF; IFNB; IFN-beta
Species	Human
Gene ID (Human)	3456
Gene ID	3456
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-187 of human IFN-beta (NP_002167.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P01574
Protein Weight	22 kDa
Shipping	Ice bag