

Anti-Ifn-Gamma Monoclonal Antibody For Wb, Fc, Elisa - P01579

Item Number: CM0006013



Introduction

Rabbit monoclonal antibody against human IFN-gamma, validated for WB, FC, ELISA. Suitable for studying immune interferon function in antimicrobial and antitumor responses. - P01579

[Learn More](#)

Parameter	Specification
Product Name	IFN-gamma Rabbit mAb
Remarks/Alias	IFG; IFI; IMD69
Species	Human
Gene ID (Human)	3458
Gene ID	3458
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-166 of human IFN-gamma (NP_000610.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, FC (intra), ELISA
Cross Reactivity	Human
SWISS	P01579
Protein Weight	19kDa
Shipping	Ice bag