

Anti-Ifn Gamma Polyclonal Antibody For Wb, Ihc-P, Elisa - P01579

Item Number: CM0013707



Introduction

Rabbit polyclonal antibody against human IFN gamma (Interferon gamma), suitable for WB, IHC-P, ELISA. Cross-reacts with human, mouse, and rat. Protein weight ~19kDa. For cytokine research and immunology studies.

[Learn More](#)

Parameter	Value
Product Name	IFN gamma Rabbit pAb
Remarks/Alias	IFG; IFI; IMD69; IFN gamma
Species	Human
Gene ID (Human)	3458
Gene ID	3458
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 24-161 of IFN gamma (NP_997468.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01579
Protein Weight	19kDa
Shipping	Ice bag