

Anti-Ifn Gamma Rabbit Monoclonal Antibody For Wb, Elisa - P01579

Item Number: CM0007403



Introduction

Rabbit monoclonal antibody targeting human IFN gamma, suitable for Western blot and ELISA. Recombinant protein immunogen; cross-reacts with human. Useful in cytokine and immune response research.

[Learn More](#)

Parameter	Value
Product Name	IFN gamma Rabbit mAb
Remarks/Alias	IFG; IFI; IMD69; IFN gamma
Species	Human
GeneID (Human)	3458
GeneID	3458
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-161 of human IFN gamma (NP_000610.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P01579
Protein Weight	19kDa
Shipping	Ice bag