

Anti-Ifn-Gamma Rabbit Monoclonal Antibody For Wb, Ip, If/Icc, Ihc-P, Elisa - P01579

Item Number: CM0011050



Introduction

Anti-IFN-gamma rabbit monoclonal antibody validated for WB, IP, IF/ICC, IHC-P, ELISA. Detects human IFN-gamma, a key cytokine in antimicrobial, antiviral, and antitumor immune responses.

[Learn More](#)

Parameter	Value
Product Name	IFN-gamma Rabbit mAb
Remarks/Alias	IFG; IFI; IMD69
Species	Human
GeneID (human)	3458
GeneID	3458
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-166 of human IFN-gamma (NP_000610.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P01579
Protein Weight	19kDa
Shipping	Ice bag