

Anti-Ifn Gamma Rabbit Polyclonal Antibody For Wb, Elisa - P01579

Item Number: CM0017179



Introduction

Offered for WB and ELISA, this rabbit polyclonal antibody targets human IFN gamma, a key cytokine in antimicrobial and antitumor immunity, supporting diagnostic and research applications.

[Learn More](#)

Parameter	Value
Product Name	IFN gamma Rabbit pAb
Remarks/Alias	IFG; IFI; IMD69; IFN gamma
Species	Human
GeneID	3458 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-166 of human IFN gamma (NP_997468.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P01579
Protein Weight	19kDa
Shipping	Ice bag