

Anti-Ifna10 Rabbit Polyclonal Antibody For Wb, Elisa - P01566

Item Number: CM0024125

Introduction

High-quality rabbit polyclonal antibody against human interferon alpha-10 (IFNA10). Validated for WB and ELISA, with cross-reactivity to rat. Ideal for antiviral research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	IFNA10 Rabbit pAb
Remarks/Alias	IFN-alphaC; IFNA10
Species	Human
GeneID (human)	3446
GeneID	3446
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-100 of human IFNA10 (NP_002162.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P01566
Protein Weight	22kDa
Shipping	Ice bag