

Anti-Ifngr1 Monoclonal Antibody For Wb, Elisa - P15260

Item Number: CM0002110



Introduction

Recombinant rabbit monoclonal antibody targeting human IFNGR1 (CD119), validated for Western blot and ELISA. Cross-reacts with mouse and rat. Ideal for studying interferon gamma signaling and immune responses.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] IFNGR1 Rabbit mAb
Remarks/Alias	CD119; IFNGR; IMD27A; IMD27B; R1
Species	Human
Gene ID (Human)	3459
Gene ID	3459
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 270-489 of human IFNGR1 (NP_000407.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15260
Protein Weight	54kDa
Shipping	Ice bag