

Anti-Ifngr1 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P15260

Item Number: CM0002065



Introduction

Rabbit monoclonal antibody targeting human IFNGR1 (CD119). Validated in WB, IF/ICC, ELISA with cross-reactivity to mouse and rat. Ideal for interferon gamma signaling and immune response studies.

[Learn More](#)

Parameter	Value
Product Name	IFNGR1 Rabbit mAb
Remarks/Alias	CD119, IFNGR, IMD27A, IMD27B, IFNGR1
Species	Human
Gene ID (Human)	3459
Gene ID	3459
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 390-489 of human IFNGR1 (P15260)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15260
Protein Weight	54kDa
Shipping	Ice bag