

Anti-Ifngr1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P15260

Item Number: CM0027274



Introduction

Rabbit polyclonal antibody against IFNGR1, validated for WB, IHC-P, ELISA in human, mouse, and rat samples. Facilitates research into interferon gamma signaling and immune response pathways.

[Learn More](#)

Parameter	Value
Product Name	IFNGR1 Rabbit pAb
Remarks/Alias	CD119; IFNGR; IMD27A; IMD27B; IFNGR1
Species	Human
GeneID (Human)	3459
GeneID	3459
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 270-489 of human IFNGR1 (NP_000407.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15260
Protein Weight	54kDa
Shipping	Ice bag